

Dear Sir/Madam

Thank you for your cooperation during the installation of new m2000 PCR system. We are shifting our assay to more sensitive technology which will decrease the lower limit of detection up to less than 12 IU/ml, for this we now required 2 gel clotted vials (serum should be 4ml) for per PCR Qualitative, Quantitative and Genotype . The sample less than the required serum volume will be rejected as QNS.

PCR department have issued the new schedule for T SPOT TB test for Lahore and Out of City.

T.SPOT TB Batch will be done only Monday, Wednesday and Friday.

All Centers which are Out of City Booking the T.SPOT TB Samples only on Sunday, Tuesday and Thursday and 100% make sure these samples should be received in head office on Batch Day before 05:00pm.

All Centers which are in Lahore City Booking the T.SPOT TB Samples only on Sunday, Monday, Tuesday, Wednesday, Thursday, Friday and 100% make sure these samples should be received in head office on Batch Day before 05:00pm.

Request to all centers mention sample collection date and time on each vial of T.SPOT TB test and ensure T.SPOT TB samples transport at Room temperature 24C- 28C and also ensure that sample must reach in the Section within 30hrs of sample collection time.

Kindly ensure that all the centers and franchise send required sample volume to avoid any inconvenience. Your Cooperation in this regard would be highly appreciated.

All regional manager will make sure to inform this message to their CLL Centers and Franchise Centers which are in your regions.

Thank you and Best Regards,

Abdul Hamid

Manager Operations

CC:

Khalida Asif
Manager Dept. of Molecular Biology
Chughtais Lahore Lab
Cell:+92 345 8445081
Ph: +92 42 35711711
khalida.asif@ccl.edu.pk